

HGVT

Org. No. A003523F

PARAFFINALIA NEWSLETTER

VOLUME 27, NUMBER 4
December 2022

The HGVT aims to provide a dynamic continuing education program in which all persons with an interest in Histology and Histotechnology are freely invited to participate.

CONTENTS:

President's Report

Reflecting on Sue Sturrock's Super Career

Under the Microscope with Pigie De Castro

Immunofluorescence for Autoimmune Enteropathy
by Bronwyn Christiansen

IHC Antibody Spotlight – CD138

Future Events 2023

Committee Page

The members of the Histology Group of Victoria and Tasmania 2020 are:

Name	Institution
Kerrie Scott-Dowell	Dorevitch Pathology/Leica Biosystems
Adrian Warmington	Ballarat OSM
Mark Bromley	Sullivan Nicolaides Pathology
Elizabeth Baranyai	Cabrini Health
Alison Boyd	Northern Health
Kellie Vukovic	Melbourne Pathology
Samantha Arandelovic	Mater Hospital Brisbane
Snejana Ursache	Alfred Hospital
Imogen Campbell	Alfred Hospital
Alistair Townsend	Royal Hobart Hospital
Cristine Gorringer	Royal Hobart Hospital
Meghan Leo	Histolab
Yashi Xie	St. Vincent's Hospital
Bronwyn Christiansen	Royal Children's Hospital
Hazel Chambers	Royal Children's Hospital
Bindi Bates	Leica Biosystems
Alex Johnston	Leica Biosystems

Contact us at

<http://www.hgvt.org.au/contact-us>

Membership enquiries: membership@hgvt.org.au

General enquiries: secretary@hgvt.org.au

Sponsorship/Advertising enquiries: trade@hgvt.org.au

Newsletter enquiries: editor@hgvt.org.au

Advertising for the next edition of Paraffinalia closes: **25th February 2023**

Paraffinalia Rates:

A4 Electronically Submitted \$400 (no GST) Flat Rate Per page
\$1200 yearly package

Positions Vacant- No Logo up to 75 words
FREE to email membership

FREE list on Website

Used Equipment- 50 words – no logos/no pictures

FREE

A4 with logo \$200 (no GST)

Articles & Reports:

Author enquiries and readers wishing to contribute articles or reports can contact the Editor - editor@hgvt.org.au

Please email articles (preferably Microsoft Word format) for inclusion in the next edition to editor@hgvt.org.au All items submitted for publication will then become the sole property of the Histology Group of Victoria Inc.

Disclaimer:

Any opinions expressed in this publication are solely those of the contributing author and are not necessarily reflective of the Histology Group of Victoria Incorporated or the editor.

NOTE: No responsibility is assumed by the Histology Group of Victoria Incorporated for any injury and/or damage to persons or property as a matter of products liability, negligence or otherwise, or from any use or operation of any methods, products, instructions or ideas contained in the material herein. It is the user's responsibility to ensure that all procedures are carried out according to appropriate Health and Safety requirements.

Copyright of this newsletter "Paraffinalia" is held by the Histology Group of Victoria Incorporated. No material may be reproduced in part or in whole without written consent from the copyright holders. All rights reserved.

President's Report – Behind the Bench

As the end to another Covid impacted year rapidly approaches, we can reflect on 2022. If you are working in a clinical histology laboratory setting, you will be buried beneath wax block trimmings and request forms. On the home front, the lawns are growing faster than you can mow, and anytime outside, you will be battling the mosquitoes. Hopefully, the rain will ease enough for us to have a real Summer and a well-earned break.

Since the last newsletter, we have had our AGM and have formed a new committee.

I would like to thank the outgoing committee for their contributions and welcome Cristine, Bron and Hazel as well as welcome back Alistair. We will begin face-to-face meetings in February at Peter Mac and will work to a State Conference later in the year.

On a personal level, I have mixed emotions around the retirement of the very special, Histo-Guru, Immuno- Queen extraordinaire, Ms Sue Sturrock. She has always been the strongest advocate for the Histology Group and further education. I am happy for her next adventure but will miss her advice and committee friendship. Good Luck Sue your contribution to the industry has been enormous.

I wish everyone the very best of the festive season. Stay safe and look forward to seeing you all in 2023.



Kerrie Scott
(Leica/ Dorevitch
Pathology)
HGVT President

ONCORETM PRO

Fully Automated Slide Stainer

RELIABLE. AUTOMATED. COMPACT.

The ONCORETM PRO Automated Slide Staining System now performs fully automated ISH so you can run your FISH probes and CISH probes with full automation and the convenience of a compact bench-top system.



The ONCORETM PRO Automated Slide Staining System is a compact and convenient bench-top instrument that is capable of performing immunohistochemistry (IHC) and *in-situ* hybridization (ISH) procedures on formalin-fixed paraffin-embedded (FFPE) tissues. The on-board capabilities include slide baking, deparaffinization, antigen retrieval, and antibody or probe detection for IHC, Multiplex IHC and ISH applications.

Contact us for more information or a quote.

www.metagene.com.au - 1800 788 498 - admin@metagene.com.au

Everything you need for a Complete Digital Pathology System

Multi Application

Comprehensive solution for a wide range of sample types, including primary H&E, IHC antibodies and FISH probes

Diagnostic Confidence

Accurate and validated computer-assisted analysis provides higher confidence in patient assessment

Workflow Efficiency

Modern computerized workflows increase lab productivity and improve FTE savings

A cutting-edge diagnostics solution, bridging the gap between Brightfield pathology and FISH. Combining Whole Slide Imaging, digital FISH analysis and revolutionary digital tissue matching of FISH with H&E/IHC samples, PathFusion provides accurate and validated analysis for higher diagnostic confidence.





PathFusion

Full Pathology Imaging Suite

H&E | IHC | FISH

POWERED BY
GENASIS



SCAN ME



Contact us for more information or a quote:

www.metagene.com.au

1800 788 498

admin@metagene.com.au

TGA APPROVED

REFLECTING ON SUE STURROCK'S

SUPER CAREER

Sue began her career in Science in a cardiac pharmacological research lab at the Austin Hospital in 1982 after being accepted into RMIT to study Medical Laboratory Science part time. As a lab technician, she prepared glassware, apparatus and solutions to study the effects of drugs on ischaemic heart disease. This might seem to be quite dry, but the opposite was the case. Comradery in the lab was fantastic and saw weekly Friday pub lunches, skiing weekends and water fights. Dungeons and Dragons was played using printer comms to and from the Uni Melb server (one of the first internet connections in Australia!) after hours, of course! While the group next door was working on the Human Genome Project, Sue watched along, observing the difficulty in development and achievement of cell culture techniques first-hand. Distilled water didn't just come out of a tap but had to be produced drip by drip from a hand-built still. She is now on one of the drugs she studied.

Knowledge attained:

- ❖ Measure cylinder volume at the bottom of the meniscus
- ❖ Rats make great friends
- ❖ How to bond with librarians over photocopiers and Index Medicus
- ❖ How to stop a lift between floors by jumping-up-and-down
- ❖ How to catch a man.

When the pay rates crossed over and a trainee wage on the Medical Scientists award was higher than the research rates, Sue obtained a position in the Histology lab at the Austin.

The decision to work in Anatomical Pathology was based purely on which lab didn't have to do night shift. The lab was also relatively new with distilled water on tap and a dishwasher (major attractions). Pathologists were aloof and distant creatures who overlapped with lab staff in the cut-up room only. There was a separate doctor's dining room. As a trainee, Sue was placed in cut-up to enter specimens in a huge log book and pass hand-written cassettes to pathologists who didn't want to be there. For sport, they would ask the newbie to hand them a 'feather' and watch as she frantically opened drawers and searched for a bird's feather. This was a strange and scary place – and then, there was - the mortuary. Releasing bodies and fetching items during autopsies was an eye opener. When eventually allowed into the histology lab, Sue discovered a witty and intelligent group of scientists, led by Tristan Roberts, who made her feel welcome and taught her everything they knew with enthusiasm and care. The wall of snow domes and postcards from all over the world and interesting personal items on display provoked educational discussions.



When the pay rates crossed over and a trainee wage on the Medical Scientists award was higher than the research rates, Sue obtained a position in the Histology lab at the Austin. The decision to work in Anatomical Pathology was based purely on which lab didn't have to do night shift. The lab was also relatively new with distilled water on tap and a dishwasher (major attractions). Pathologists were aloof and distant creatures who overlapped with lab staff in the cut-up room only. There was a separate doctor's dining room. As a trainee, Sue was placed in cut-up to enter specimens in a huge log book and pass hand-written cassettes to pathologists who didn't want to be there. For sport, they would ask the newbie to hand them a 'feather' and watch as she frantically opened drawers and searched for a bird's feather. This was a strange and scary place – and then, there was - the mortuary. Releasing bodies and fetching items during autopsies was an eye opener. When eventually allowed into the histology lab, Sue discovered a witty and intelligent group of scientists, led by Tristan Roberts, who made her feel welcome and taught her everything they knew with enthusiasm and care. The wall of snow domes and postcards from all over the world and interesting personal items on display provoked educational discussions.

The department tearoom eventually became a mecca for pathologists, technical and scientific staff after a brave young pathologist, Bill Murray, breached the divide. Even surgeons were seen in later years enjoying the ambiance and the magnificent cakes.

Nobody was allowed to use anyone else's steel knife and the quality of the knife was related to your status and experience. If you cut a bone or a calcified specimen you had to sharpen your knife again. You sharpened these dangerous heavy pieces of steel on a Tem Tool – a flat circular plate which whizzed around at high speed onto which you lowered the knife. It ground the nicks off the knife's edge or, if not screwed on properly, was flung across the room. Stropping the edge to remove rough particles was the final stage. Needless to say, knives were treated like gold and were recognised as the most dangerous and delicate item in the lab. If you dropped a knife you jumped out of the way. It could cut your foot off! When disposable knives were introduced it was the dawn of a new age.



Sue witnessed the development and introduction of immunohistochemistry at the Austin in the mid 1980's. From handstaining in Perspex racks to Sequenza units and finally, automated stainers, Sue saw the span of antibodies increase from CEA and CAM5.2 to 160 different targets by the late 2000's. As the AIDS epidemic took hold, she developed histochemical and IHC stains for the terrible range of infections these poor people grew. A tragedy for the individuals, HIV became a smorgasbord of scientific

interest for Histology and Microbiology laboratories, resulting in presentation of findings at the AIMS conference in Melbourne. The Austin was developing adult liver transplants at the same time and Sue participated in the assessment of donor livers prior to transplant via early morning Oil Red O's on frozen section liver samples (so much for no night work!). The transplant physicians were a lovely bunch, but you had to drop everything when they called. Worst ever was when I had just finished the supermarket shopping and while at the checkout, had to leave it all and run! Renal biopsy attendance and management saw Sue develop skills in glomerulus spotting, resin processing and glass knife breaking. When this was deemed too slow, paraffin processing saw her in a walk-in fridge cutting them and spending a significant amount of time in the electron microscopy suite where similar cool temperatures were perfect for wine storage. Wine education was gained here and with a group of like-minded fellow scientists monthly at the Austin and the newly formed HGV and satellite IHC groups.

The HGV was a career highlight meeting many people and sharing scientific experiences at committee meetings and conferences organised. Newsletter collation nights were hard but fun, done by hand and sent via snail mail. Histologists are great people always ready to help anyway they can.

In the 1990's, histology was mostly about diagnosing cancer. It was a celebration when a patient was diagnosed with a benign condition. When the Ludwig Institute formed a collaboration with the Austin to identify targets on tissue for antibody therapy for solid tumours, Sue and colleagues created a Tissue Bank of normal and tumour tissues and used IHC to screen chimeric (mouse-human) antibodies for their suitability as carriers of targeted therapies. This was the beginning of immunotherapy and although only one suitable colonic carcinoma antibody was identified, it showed the promise of the individual targets and therapies IHC could identify for cancer.



Events making an impression:

- ❖ Watching associates develop PCR methodology
- ❖ Roof tiles falling in as a result of water pipes breaking in floor above, simultaneously with power out as senior scientist held an umbrella over paper files
- ❖ Celestine blue boiling over on hotplate (as soon as you turned your back)
- ❖ Pan flusher leaking through ceiling at special stains
- ❖ Dirty Dancing
- ❖ Wax crack impression of different senior scientist
- ❖ Pink Schiff's hands whenever a social event was planned
- ❖ Chloroform-affected ski boot fitting
- ❖ Chasing bowel down plughole-slippery suckers
- ❖ Gourmet Christmas lunches by invitation only and goodie box.

With 2 children produced and a Northcote weatherboard renovated it was time to look for a career change for Sue at Peter MacCallum in 2007, re-uniting with Bill Murray, now Head of Anatomical Pathology and Louise Keelan in Cytology. Going from a large pathology department to a small one, re-established only 5 years earlier, was certainly a change of pace, but provided significant challenges. The lab was tiny, the size of a lounge room, but fitted microtomy, special stains, cytology and IHC into one area. Everyone's personal space was reduced, but the window to the Fitzroy Gardens provided welcome perspective. CISH, to detect Her2 breast cancer, was run by hand and was the first of many novel assays to be established. This was a major money spinner allowing the Molecular laboratories to be grown, a renovation of the cut-up area and processor and stainer replacement.

As the lab grew, Sue was lucky to engage with wonderful, technically excellent Scientists complemented by many RMIT professional practice students. FISH was established to detect ALK, a translocation marker of eligibility for a very successful drug therapy, stopping lung cancer in its tracks for people who had never smoked. This was very personal for Sue and the other scientists in the group, as it was the first time histologists could feel part of a cancer cure, rather than just providing bad news. With changes to NPAAC guidelines in 2013, a complex cutup training programme for scientists was commenced with training by Canadian pathologists' assistant, Courtney Savill. This created another career direction for scientists and released pathologists to do what they do best – microscopic diagnosis.

The greatest challenge of Sue's career was moving the Peter Mac lab to its current location in Parkville in 2016. Supported by the hospital's long lead into the change process and the great team in the lab, the move was very successful and produced a much bigger space for both histology and cut-up. The best planning and architects in the world could not foresee the car dashes between both sites to save bench supports for the new, worryingly sagging microtome benches or the failure to turn on the deionised water system. The new site produced further connections with pharmaceutical companies and more rapid breakthroughs in cancer therapy, with IHC again providing screening for confirmation by molecular assays and FISH for melanoma, lung, ovarian, colon and other cancers.

Peter Mac Highlights:

- ❖ Daily commute to the Commune for morning tea
- ❖ Managing deionized water to be directly tapped to the lab (career-long concern)
- ❖ Free entry to the last session of the Boxing Day tests at MCG
- ❖ Long Island teas at roof top bars in Summer
- ❖ Meeting Carlton footballers and Penny Wong at Parliament
- ❖ Watching a cascade of block boxes fall off collapsing shelving onto a pathologist's desk and the months of sorting afterwards
- ❖ NSH conference in Austin, Texas
- ❖ Central Melbourne dinners with wine and guilt free train travel home
- ❖ Chased by a possum in Fitzroy Gardens
- ❖ Same day lobster.

Looking for a change of pace after the tumultuous movement of staff and location at Peter Mac, Sue began at Melbourne Pathology in the histology lab in 2018. Having mostly worked in the public sector, the move to private pathology provided a welcome change despite the huge volume of work. In a return to her roots, the pre-established economy of procedure and provision of top-notch equipment produced a simplicity allowing the histology to be the main focus. As trainer and problem-solver, Sue was re-invigorated and relished the opportunity to pass on her knowledge to students and learn from pathologists and scientists in the group. Adapting to the Covid experience showed the calm, rapid capability for versatility at Melbourne Pathology. All employees retained their positions and many more were employed to cope with the racks and racks of Covid tests. Everyone from the CEO down helped batch in a fantastic team effort.

Greatest advances seen:

- ❖ Disposable blades
- ❖ Barcode driven tracking for cassette and slide labelling
- ❖ Predictive IHC markers for therapeutic cancer treatment

Sue will remain grateful to the many who have supported her through her pathology career.





Artisan Link Pro

The consistent, safe and easy choice for special stains

The Artisan Link Pro Special Staining solution is designed to provide pathology labs with consistent staining results, optimal user safety and significantly improved efficiency. Automate your special stains with Artisan Link Pro and experience an easy and uncomplicated process for consistent results, slide after slide.

- Up to 14 protocols stained simultaneously
- 28 high-quality Ready-to-Use standardized stain kits
- 30 optimized, validated protocols that are editable for easy, yet flexible, results

www.agilent.com

© Agilent Technologies Denmark ApS, 2018



TekEquipment

Equipment for a safer, cleaner & more efficient laboratory.

Down-Draft GrossPath GP Workstations

- Integrated active carbon air recirculation filter (+ HEPA option)
- No ducted ventilation system needed
- Height adjustable
- Simple installation
- Full range of accessories
- Widths available: 125cm, 150cm, 175cm, 200cm
- Cirdan PathLite Digital Camera system & PC option



Cut Grinder Diamond Pathology Saws

- Water cooling of blade eliminates aerosol material
- Safe to touch diamond coated cutting band
- Quick tensioning & variable speed band
- Perfect cut surface
- Full stainless steel construction
- 22 or 31cm cutting height models available



ErgoTek

Ask us about our specially designed ergonomic benchtop!

With 50cm of height adjustment, a 400kg capacity, a range of colour options available & rock solid stability with microtome use in mind, it's a must have for your laboratory.

Oscillating Autopsy Saws

- Variable speed oscillation & gradual start up to protect mechanical components
- Easy to clean & simple to change saw blades
- Slim ergonomic construction designed for long periods of use
- Quick charging battery powered options available



Ventilated Tissue Storage Cabinets

- Ducted or integrated fan & filter
- Stainless steel or glass doors
- Stainless steel or coated laminate construction
- Height adjustable shelving
- Roll out trolley option



Contact Us today to learn more!

1300 368 138 | tekequipment.com.au | info@tekequipment.com.au



UNDER THE MICROSCOPE

with Pigie De Castro

What is your current job?

I currently work as a Medical Scientist at The Royal Children's Hospital

How long have you been working in Histology?

I did my work-placement in 2008, then the year after that I was employed as a trainee scientist while I finish my last year of university part time. Since then I have always worked in Histo. So roughly 13 years.

What piece of advice would you give to people starting in Histology?

Learn as much as you can and don't be afraid to ask questions. Also learn/watch how others work and pick up the skills/techniques that you think works well for you.



What's a goal you have for this year?

Start and finish refurbishing the wooden furniture I have hoarded! Our garage looks like an old furniture storage room. Lol

What music is on your playlist at the moment?

I currently have Ed Sheeran radio playlist on my Spotify account.

Other than histology, what is another scientific topic that interests you?

I am into veggie gardening at the moment and I find soil biology very interesting. Learning about soil food web and how it helps a plant grow is mind boggling.

What movie or TV show do you think everyone should watch?

Jason Bourne movies. I would always sit and watch them from start to finish every-time I see it.

If you could live anywhere in the world for a year, where would you pick?

I would choose Norway because it is a beautiful country. I would visit all the Fjords there.

If you could choose two famous people, past or present, to have dinner with, who would it be?

Sir David Attenborough and Alex Hannold. I would love to hear the stories about all the places they have been and adventures that they had.

VEGA Laser Cassette Marker




The Vega Laser Cassette Marker offers permanent indelible marking of tissue cassettes with text, linear barcodes, and two-dimensional barcodes.




VEGA FEATURES:

- **Unique phase technology permanently marks cassette surfaces**
Cassettes are printed using an innovative low-powered YAG laser. Text and barcode images are high resolution (600 dpi)
- **Small footprint, fast print speed, high resolution**
12 cassettes per minute (one every 5 seconds)
- **Six cassette magazines**
Up to six different cassette colours simultaneously
- **Full LIS/LIMS system integration**
Supports linear and two-dimensional barcode and true type fonts
- **Reduced operational expenses**
Eliminates the need for ink cartridges, ribbon and UV fixing lights



DM 100 Characters



DM 50 Characters

Trusted by Australian laboratories with a team of dedicated service & support staff in major capital cities



Distributed by Abacus dx - Contact Mick for a quote


Mick van Werkhoven

Territory Manager- VIC & TAS

0499 400 601 | m.vanwerkhoven@abacusdx.com | www.abacusdx.com

abacus dx

Pyramid Innovation's Smart Printer Range



Scan & Print with ease

Track & Trace your lab workflow in real time



*Twin Hopper
Slide Printer*



*Six Hopper
Cassette Printer*



Print
standard
slides,
IHC adhesion,
specials like ThinPrep
or control slides
in the same print job
in the right order!

Pedigree

Built by the team that invented
the industry standard Microwriters.

Performance

Fast - 4 second thermal labelling for
cassettes & microscope slides.

Future-proof

Not 'just a printer' - Full Track & Trace
suite to digitise your workflow with
lifetime warranty options available.

Thin & LEAN

Perfect for labs starved of bench space.

Touchscreen

Quality of life features at your fingertips.

Tracking

Complete chain of custody handling.
Rapidly search and report on your lab
in real time - live as you work.

Scanner

Scans all forms of barcodes including
2D, linear and QR code.

LIS / On Demand

Compatible with all LIS / LIMS as well as
offering true 'On Demand' LEAN.

Stellar Service

Included preventative maintenance &
lifetime product support with additional
training, customisation and live support
when you and your staff need it.
Find out why our clients love us.

Scan to
find out more &
arrange a demonstration



IMMUNOFLUORESCENCE FOR AUTOIMMUNE ENTEROPATHY

by Bronwyn Christiansen


Autoimmune enteropathy is caused by an autoimmune attack. It usually presents in the first six months of life and principally involves the small intestine, although the stomach and large intestine may also be involved. Autoimmune enteropathy is a rare cause of intractable diarrhea associated with circulating gut autoantibodies and a predisposition to autoimmunity.

Symptoms of this disorder may include diarrhea / loose watery stools all the time, poor weight gain and weight loss, decreased urine output, frequent infections, occasional blood in the stool and skin rash. Autoimmune enteropathy may be linked with other disorders like diabetes and kidney disease.


For children with suspected anti enterocyte/anti-goblet cell enteropathy, an endoscopic procedure is performed to gain a diagnosis.

Anti-enterocyte brush border antibodies are demonstrable by indirect immunofluorescence.

An indirect immunofluorescence technique can be used to identify circulating auto antibodies to mucosal structures in a patient's serum. The patient's serum is incubated with frozen sections of the substrate (a normal donor tissue – commonly frozen duodenum) then the bound auto antibodies are detected with the use of anti-human immunoglobulins.



An example of anti-goblet cell antibody binding after incubation with patient serum followed by IgA/G FITC.



An example of anti-enterocyte antibody binding after incubation with patient serum followed by IgA/G FITC. Note the brush border staining and negative nuclei.

Treatment

The treatment of this problem is usually by medicines that suppress the immune system. Children may also need a special diet. Surgery is usually not needed. Sometimes these children will need to get IV (intravenous) nutrition. This gives them the nutrition their bodies need to grow and heal the intestines.

References

Adult Autoimmune Enteropathy: Mayo Clinic Rochester Experience, Akram S et al. Clinical Gastroenterology and Hepatology. 2007; 5: 1282-1290.



Experience the future of digital pathology

Prepare your laboratory for the future with digital pathology solutions from 3DHISTECH™ and Epredia™

- Whole-slide scanners that accurately capture minute details at incredible speeds
- Software systems that enable collaboration across laboratories – and industries
- Microarrays enabling high-efficiency processing and high-density storage



Pannoramic 1000

When speed is your need. Free up your lab's resources with the new Pannoramic 1000 – the high-throughput whole-slide scanning solution.



Pannoramic 250 FLASH III

3DHISTECH Pannoramic 250 Flash III, an all-in-one solution or digital pathology research and storage. Enjoy increased speed and efficiency.



Pannoramic SCAN II

Save time in routine pathology; both brightfield and fluorescent scanning solution in the same machine.



Pannoramic MIDI II

A versatile, low-volume digital pathology solution for smaller labs.



Pannoramic DESK II

An excellent choice for tele-consultation and remote section scanning.



TMA Grand Master

High-capacity workflow; high-speed microarray; simultaneous loading, imaging, drilling & punching



Bio-Strategy PTY Limited

T 1800 00 84 53 | sales.au@bio-strategy.com

www.bio-strategy.com | shop.bio-strategy.com



bio-strategy
delivering technology

BOND-PRIME

TRANSCEND TRADITIONAL WORKFLOWS WITH **UNIVERSAL ACCESS**

Seamlessly adapt to incoming workflows – Batch, Continuous, Single Slide.

Only the new BOND-PRIME IHC system delivers all elements of Universal Access.



LOAD ANY SLIDE

Load any staining protocol, chromogen, or technology.



IN ANY COMBINATION



Load single slides, complete cases or full batches.



ACTIVE REAGENT CONTROL (ARC) MODULES

Individual staining chambers allow you to run any staining protocol at any time.



WITH ANY REAGENT



Add markers or detection at any time with 70 reagent positions.

AT ANY TIME



The Preload & Unload Drawers let you load & unload slides on your schedule, not the instrument's.



BOND-PRIME ADAPTS TO ANY WORKFLOW:



Batch



Continuous



Single slide

ACHIEVE 100% INSTRUMENT UTILIZATION AND RELEASE YOUR FULL DIAGNOSTIC PRODUCTIVITY

Advancing Cancer Diagnostics
Improving Lives

FOR IN VITRO DIAGNOSTIC USE


Leica
BIO SYSTEMS

Copyright © 2022 Leica Biosystems, a division of Leica Microsystems Inc. All Rights Reserved. LEICA and the Leica logo are registered trademarks of Leica Microsystems IR GmbH. BOND-PRIME is a trademark of the Leica Biosystems group of companies in the USA and optionally in other countries. Other logos, product and/or company names might be trademarks of their respective owners.

220275 Rev A • 04/2022

IHC ANTIBODY SPOTLIGHT: CD138

One of the most abundant surface proteins on the surface of myeloma and plasma cells is syndecan-1 or CD138. (Akhmetzyanova et al., 2019) The CD138 molecule is a transmembrane heparan sulphate glycoprotein expressed at distinct stages of differentiation in normal lymphoid cells such as pre-B cells, immature B cells and Ig-producing plasma cells as well as being expressed in stratified and simple epithelia. The loss of CD138 expression from atypical cells is reported to be an early event during cervical carcinogenesis whereas CD138 antigen expression shows a close association with preserved epithelial morphology and differentiation. However, CD138 is commonly used to identify and quantitate plasma cells. In normal tissues, CD138 is expressed on plasma cells but also in various epithelial cell types. CD138 is also expressed in various cancers. In several tumor types CD138 expression levels were described to be prognostically relevant. CD138 expression in cancer is of potential clinical interest. Antibody-based drugs targeting CD138 in plasmacytomas are being evaluated in clinical trials. In preclinical studies, anti-CD138 antibodies were also effective against triple negative breast cancer and melanoma cells.



[LEFT] Breast ductal carcinoma: Breast cancer no special type, NST) with strong membranous CD138 staining of tumor cells and intensive CD138 positivity of the tumor stroma

[RIGHT] Plasmacytoma: Note the intense membrane staining of tumor cells

References

Akhmetzyanova, I. et al. (2019) Dynamic CD138 surface expression regulates switch between myeloma growth and dissemination, Nature News. Nature Publishing Group. Available at: <https://www.nature.com/articles/s41375-019-0519-4> (Accessed: December 4, 2022).

CD138 (syndecan 1) (2022) CD138 (Syndecan 1) - IHC Primary Antibodies. Available at: <https://shop.leicabiosystems.com/en-au/ihc-ish/ihc-primary-antibodies/pid-cd138-syndecan-1> (Accessed: December 4, 2022).

Kind S. et al. (2019) A shift from membranous and stromal syndecan-1 (CD138) expression to cytoplasmic CD138 expression is associated with poor prognosis in breast cancer. Mol. Carcinogenesis, 58(12):2306-2315

Kind, S. et al. (2019) Prevalence of Syndecan-1 (CD138) Expression in different kinds of human tumors and normal tissues. Disease Markers, Volume 2019, Article ID 4928315



Hello Histo group! The HGVT is interested in promoting Histology into the public eye more (to show off of course!) and want to see pictures from our community to develop an Instagram (or Histo-gram, for the pun-minded). Above is some artistic inspo of the already incredible potential our discipline offers, however, anything related to Histology, your lab or your team is more than welcome!

Send all pictures to: editor@hgvt.org.au

Best regards!

Links to above art:

<https://www.nsh.org/blogs/christine-carreira/2021/03/26/art-of-the-stain-2021-contest-winner>

<https://www.pinterest.com.au/pin/20055160831353005/>

<https://illnessnarratives.files.wordpress.com/2014/01/i-will-wear-my-heart-upon-my-sleeve.jpg>

<https://www.pinterest.com.au/pin/395331673514803840/>

BenchMark ULTRA PLUS

Experience the power of PLUS



Optimised workflow ▪ Flexible solutions ▪ Quality is at the CORE

- **New proprietary waste management solution** degrades DAB to less than 5 ppb, meeting most stringent municipality requirements for nonhazardous waste.
- **VSS 14.x software** allows for easier navigation of inventory management and ordering through a new filter as well as lot-to-lot validation on the same instrument.
- **New integrated touchscreen** reduces bottlenecks by allowing many functions to be done at the instrument.
- **New slide drawer design** mitigates fluid ingress and egress resulting in continuous and highly reliable drawer performance.
- **New retractable work surface** gives technicians additional space to work while at the instrument.

**THIS PRODUCT IS NOT AVAILABLE FOR PURCHASE BY THE GENERAL PUBLIC.
ALWAYS READ THE LABEL AND FOLLOW THE DIRECTIONS FOR USE.**

© 2022 Roche Diagnostics. All rights reserved.

BENCHMARK is a trademark of Roche.

Roche Diagnostics Australia Pty Limited 2 Julius Avenue North Ryde NSW 2113
Ph: +61 2 9860 2222 ABN 29 003 001 205
MC-AU-02065



EUKITT® Mounting Medium

Protecting what is most important to you

Q: Do you find it difficult to source a high quality mountant?

EUKITT®

A: EUKITT® mounting medium provides the ultimate protection to ensure precious samples are well preserved for years to come.

EUKITT® has long been regarded as the gold standard of mounting media for histology. EUKITT® is characterised by fast curing, crystal-clear optics with a favourable refractive index, good fluidity, low shrinkage, and can be used manually and in automated coverslipping equipment.

Resistant to discolouration, heat, cold, humidity, daylight and UV light, laboratories all around the world trust EUKITT® in the preparation of their specimens.



- Refractive index n [20 °C]: 1.49
- Viscosity: 250 – 450 mPa*s
- Curing time: 20 min
- No yellowing after long storage
- No brittle fracture
- No primary fluorescence under UV light

Showing the right path

Contact us on 1800 257 213 to learn more about how EUKITT® mounting medium can protect your specimens.

www.trajanscimed.com



Future Events: 2022

Org. No. A0035235F

Date: 23rd February, 2023

Scientific Meeting

Topic: Short Presentations

Venue: Zoom Meeting (streamed and recorded)

Date: 20th April, 2023

HGVT Scientific Meeting

Topic: TBA **Venue:** TBA